

RevertUP™ II Reverse Transcriptase, 200 U/µl

LOT: See product label EXPIRY DATE: See product label

ORDERING INFORMATION

CAT. NO.	SIZE	PACKAGE CONTENT
BR0400501	10000 U	50 μl RevertUP II Reverse Transcriptase 1 ml 5× Reverse Transcriptase Buffer
BR0400502	50000 U	2 × 125 μl RevertUP II Reverse Transcriptase 2 × 1 ml 5× Reverse Transcriptase Buffer

COMPONENT	COMPOSITION
RevertUP II Reverse Transcriptase	RevertUP II Reverse Transcriptase, 200 U/µl, in storage buffer containing 50% (v/v) glycerol.
5× Reverse Transcriptase Buffer	Optimized 5× Reverse Transcriptase reaction buffer.

STORAGE -20°C (until expiry date – see product label)

FEATURES

- Improved performance for synthesis of long cDNAs (≥ 19 kb)
- Excellent stability at high temperatures up to 60°C
- High sensitivity for cDNA synthesis from few copies of template

APPLICATIONS

- First-strand cDNA synthesis
- Generation of labeled cDNA
- RNA analysis by primer extension
- cDNA library construction
- RT-PCR

DESCRIPTION

biotechrabbit[™] RevertUP II Reverse Transcriptase is a proprietary MMuLV reverse transcriptase engineered by point mutations resulting in increased temperature stability up to 60°C. This guarantees top performance with templates showing a high degree of secondary structures. The enzyme has no RNase H activity and ensures efficient synthesis of ≥19 kb cDNA.

RevertUP Reverse Transcriptase is a DNA polymerase which uses RNA as a substrate and exhibits no measurable proofreading $3' \rightarrow 5'$ exonuclease function. This enzyme performs cDNA synthesis by extending a DNA primer annealed to an RNA template; it can also make copies of single-stranded DNA templates.

PROTOCOL

Prevention of cDNA synthesis reaction contamination

RNase contamination is an exceptional concern when working with RNA. RNase A, providing most threat to RNA integrity, is a highly stable contaminant of any laboratory. To prevent RNA from degradation and to minimize possibility of contamination during cDNA synthesis; follow the guidelines below:

- Use separate clean areas for preparation of the samples and the reaction mixture.
- DEPC-treat all tubes and pipette tips or use certified nuclease-free labware with aerosol filters.
- Wear fresh gloves when handling RNA and all reagents.
- Always assess the integrity of RNA prior to cDNA synthesis in denaturing agarose gel electrophoresis.
- Use RNase free water and other reagents.
- To prevent RNA from degradation, add Ribonuclease inhibitor (optional) in to the cDNA synthesis reaction (20 units for 20 µl reaction).

Typical cDNA synthesis reaction set up

- Thaw on ice and mix very well all reagents.
- Assemble and keep all reactions on ice.
- To use time and reagents effectively, always prepare master mix for multiple reactions. For a master mix volume, always calculate the number reactions that you need plus one additional.
- Combine the following in an RNase-free reaction tube:

COMPONENT	VOLUME	FINAL CONCENTRATION
dNTP Mix (10 mM each dNTP)	2 µl	1 mM (each dNTP)
RNase Inhibitor, 40 U/µl (optional)	0.5 µl	1 U/µl
Oligo (dT) ₁₂₋₁₈ (10 μM) – or	0.5 µl	0.25 μM
Hexamer Primer (25 µM) – or	1 µl	1.25 μM
Gene Specific Primer (10 µM)	0.5 µl	0.25 μΜ
5× Reverse Transcriptase Buffer	4 µl	1×
RNA Template	0.1–1 µg total RNA or 50–500 ng mRNA (polyA)	
RevertUP II Reverse Transcriptase	1 µl	10 U/µl
RNase-free water	Variable	
Total volume	20 µl	

- Mix and collect the drops by centrifuging briefly.
- When using
 - Hexamer Primer, incubate 10 minutes at 30°C followed by 42–55°C for 20–60 minutes
 - Oligo (dT) or gene-specific Primer incubate at 42–55°C for 20–60 minutes.
- Inactivate enzyme at 99°C for 5 minutes.
- Collect the drops by spinning briefly.
- Store products at -20°C or proceed to next step, like PCR or qPCR.
- Use maximum 10 µl of the cDNA synthesis reaction mix for PCR in 50 µl volume.

CERTIFICATE OF ANALYSIS

Unit Definition

One unit is the amount of enzyme activity that incorporates 1 nmole of dTTP into acid insoluble fraction in 10 minutes at 42°C when poly(A)+ RNA and oligo(dT)20 are used as template–primer.

Quality Control

Exonuclease assay

Linearized lambda/Hindll fragments are incubated with the Reverse Transcriptase in a 50 µl reaction mixture for 4 h at 37°C. No degradation of DNA was observed.

Endonuclease assay

lambda DNA is incubated with the Reverse Transcriptase in a 50 μ l reaction mixture for 4 h at 37°C. No degradation of DNA was observed.

Nick Activity

Supercoiled plasmid DNA is incubated with the Reverse Transcriptase in a 50 µl reaction mixture for 4 h at 37°C. No conversion of covalently closed circular DNA to nicked DNA was detected.

E. coli genomic DNA contamination assay

A sample of Reverse Transcriptase is analyzed with specific primers targeting the 16S rRNA gene in qPCR for the presence of contaminating *E. coli* DNA. No *E. coli* DNA was detectable.

RNase Assay

A sample of the enzyme was incubated with a RNA template. RNase activity was not observed after agarose gel electrophoresis.

Functional Assay

cDNA synthesis with specific primers, followed by quantitative PCR.

Quality confirmed by: Head of Quality Control

SAFETY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter: http://www.biotechrabbit.com/support/documentation.html.

USEFUL HINTS

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

CONTACT BIOTECHRABBIT

biotechrabbit GmbH Volmerstr. 9a info@biotechrabbit.com 12489 Berlin, support@biotechrabbit.com Germany www.biotechrabbit.com

Phone: +49 30 555 7821-10 Fax: +49 30 555 7821-99



Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.

Trademarks: biotechrabbit[™], RevertUP[™] (biotechrabbit GmbH).

valid from 13.09.2019

Certified Quality-System

mdc

ISO 13485