

25 mM dNTP Mix

LOT: See product label EXPIRY DATE: See product label

ORDERING INFORMATION

CAT. NO.	SIZE	PACKAGE CONTENT
BR0600502	1 ml	1 ml 25 mM dNTP Mix (25 mM dATP, 25 mM dCTP, 25 mM dGTP, 25 mM dTTP)
BR0600503	5 ml	5 × 1 ml 25 mM dNTP Mix (25 mM dATP, 25 mM dCTP, 25 mM dGTP, 25 mM dTTP)

COMPONENT	COMPOSITION
25 mM dNTP Mix	Aqueous solution (pH 7.0) containing 25 mM 2'-deoxyadenosine 5'-triphosphate sodium salt, 25 mM 2'-deoxycytidine 5'-triphosphate sodium salt, 25 mM 2'-deoxyguanosine 5'-triphosphate sodium salt, 25 mM 2'-deoxythymidine 5'-triphosphate sodium salt
STORAGE	-20°C (until expiry date - see product label)

FEATURES

- Exceptional quality dNTPs of >99% purity confirmed by HPLC
- Free from DNA and PCR inhibitors
- Consistent PCR results due to outstanding dNTPs stability

APPLICATIONS

- Standard or hot-start PCR
- · Long-range and high-fidelity PCR
- cDNA synthesis and RT-PCR
- qPCR
- Sequencing
- DNA labeling

DESCRIPTION

biotechrabbit[™] deoxynucleotide triphosphates are first-choice nucleotides for all PCR applications, including the most demanding, such as amplification of long targets (up to 40 kb), GC-rich templates, qPCR, cDNA synthesis, high-fidelity PCR, DNA labeling and sequencing.

Advanced production technology ensures that deoxyribonucleotide triphosphates have >99% purity and outstanding stability, ensuring excellent performance and consistent, reliable results.

For the maximum flexibility, nucleotides are available in sets and mixes of common concentrations.

PROTOCOL

Prevention of PCR contamination

When assembling the amplification reactions, care should be taken to eliminate the possibility of contamination with undesired DNA.


- Use separate clean areas for preparation of samples and reaction mixtures and for cycling.
- Wear fresh gloves. Use sterile tubes and pipette tips with aerosol filters for PCR setup.
- Use only water and reagents that are free of DNA and nucleases.
- With every PCR setup, perform a contamination control reaction that does not include template DNA.

Typical use of dNTPs


- 0.2 mM of each nucleotide is the recommended final dNTP concentration in common PCR applications.
- Every time before the use thaw and mix well each dNTP solution, especially highly concentrated ones.
- Use the ready dNTP mixes in PCR (at 0.2 mM final concentration each) according to the guidelines below:

		Volumes of dNTP mixes for the indicated PCR volume			
dNTP Mix concentration	Used for PCR as	50 µl PCR	25 µl PCR	20 µl PCR	
25 mM	125×	0.4 µl	0.2 µl	0.16 µl	

DEOXYRIBONUCLEOTIDE CHEMICAL FORMULA AND STRUCTURAL FORMULA TRIPHOSPHATE MOLECULAR WEIGHT



2'-deoxythymidine 5'triphosphate C₁₀H₁₄N₂O₁₄P₃Na₃ MW 548.1 (acid form 482.1)



25 mM dNTP Mix

CERTIFICATE OF ANALYSIS

Quality Control

Functional assay Functionality tested in PCR.

Purity test

The concentration of each lot is verified by optical density spectrometry. The purity of nucleotide triphosphates is determined by HPLC. The minimum passing specification for the triphosphate content is 99%.

Contamination test The absence of human and *E. coli* DNA is confirmed by qPCR.

Quality confirmed by: Head of Quality Control

SAFETY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter:

http://www.biotechrabbit.com/support/documentation.html.

USEFUL HINTS

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- · Most biotechrabbit products are available in custom formulations and bulk amounts.

CONTACT BIOTECHRABBIT

biotechrabbit GmbH Volmerstr. 9a info@biotechrabbit.com 12489 Berlin, support@biotechrabbit.com Germany www.biotechrabbit.com

Phone: +49 30 555 7821-10 Fax: +49 30 555 7821-99



Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.

Trademarks: biotechrabbit™ (biotechrabbit GmbH).

valid from 16.05.2019

Certified Quality-System

ISO 13485