

ORDERING INFORMATION

PRODUCT	GenUP™ Total RNA Kit		
CAT. NO.	BR0700901	BR0700902	BR0700903
SIZE	10 preps	50 preps	250 preps
COMPONENTS			
Buffer LYSIS LR	15 ml	30 ml	160 ml
Buffer WASH A (concentrate)	5 ml (add 5 ml ethanol)	15 ml (add 15 ml ethanol)	70 ml (add 70 ml ethanol)
Buffer WASH B (concentrate)	6 ml (add 24 ml ethanol)	16 ml (add 64 ml ethanol)	36 ml (add 144 ml ethanol)
Water, RNase-free (for ELUTION)	2 ml	3 × 2 ml	25 ml
Mini-Filter (blue)	10	50	5 × 50
Mini-Filter (violet)	10	50	5 × 50
Collection Tubes (2 ml)	50	5 × 50	25 × 50
Elution Tubes (1.5 ml)	10	50	5 × 50

STORAGE

Room temperature (until expiry date – see product label).
If precipitation appears, gently warm the solution to dissolve the precipitate.

FEATURES

- Fast and simple procedure
- High yields of pure RNA
- Physical removal of DNA, no DNase treatment, no toxic β -mercaptoethanol

APPLICATION

- Universal kit for total RNA isolation from various sources and different amounts of starting material

GenUP™ Total RNA Kit

DESCRIPTION

biotechrabbit™ GenUP Total RNA Kit has been specially developed for a quick and easy purification of total RNA from eukaryotic cell suspensions, tissues and biopsies, Gram-negative (e.g., *E. coli*) and Gram-positive bacteria and other sources. After few initial procedures, the RNA is bound to a filter, washed and then eluted in a separate tube. DNA is removed physically by binding to a filter without any DNase treatment or the use of toxic β -mercaptoethanol. The purified RNA is ready to be used in all demanding molecular biology applications, including cDNA synthesis, northern blot analysis and others.

SPECIFICATION

STARTING MATERIAL	Eukaryotic cells (5×10^6), tissue samples (up to 20 mg), bacterial cells (Gram-positive or Gram-negative, 1×10^9)
EXTRACTION TIME	Approximately 20–40 min
BINDING CAPACITY	100 μ g RNA
TYPICAL YIELD	Yield is highly dependent on sample type

MATERIALS SUPPLIED BY THE USER

- 70 % ethanol (molecular biology grade, non-denatured)
- 96–99.8 % ethanol (molecular biology grade, non-denatured)
- Centrifugation tubes
- Pipette tips
- *Optional:* DNase I
- For bacteria: TE Buffer: (10 mM Tris HCl, 1 mM EDTA; pH 8.0) in nuclease free water
- For bacteria: 50 mg/ml lysozyme or other bacterial lysis protein in nuclease free water

STEPS BEFORE STARTING

- Add the following volume of 96–99.8 % ethanol to each buffer bottle, close firmly, mix thoroughly and store at room temperature.

CAT. NO.	CONCENTRATE	ETHANOL	FINAL VOLUME
Buffer WASH A	BR0700901	5 ml	5 ml
	BR0700902	15 ml	15 ml
	BR0700903	70 ml	70 ml
Buffer WASH B	BR0700901	6 ml	24 ml
	BR0700902	16 ml	64 ml
	BR0700903	36 ml	144 ml

- Perform all centrifugation steps at room temperature.
- Avoid repeated freezing and thawing of frozen samples.
- Mark all vials and filters to avoid confusion when purifying multiple preps.

GUIDELINES FOR PREVENTION OF RNA DEGRADATION

Special care should be taken to minimize contamination with RNases, as RNA is extremely sensitive to degradation.

- Always wear gloves and change them frequently.
- Keep all tubes closed when possible.
- Keep samples and isolated RNA on ice.
- Reduce preparation time as much as possible.
- Use only sterile, disposable polypropylene tubes throughout the procedure (these tubes are generally RNase-free).
- Non-disposable plastic ware should be treated before use to ensure that it is RNase-free. Plastic ware should be thoroughly rinsed with 0.1 M NaOH, 1 mM EDTA followed by RNase-free water. You can also take chloroform-resistant plastic ware rinsed with chloroform to inactivate RNases.
- All glassware should be treated before use to ensure that it is RNase-free.
 - Glassware should be cleaned with detergent, thoroughly rinsed and oven baked at 240 °C for four or more hours before use. Oven baking inactivates RNases and ensures that no other nucleic acids (such as plasmid DNA) are present on the surface of the glassware.
 - Autoclaving alone will not inactivate many RNases completely. The glassware should be immersed in 0.1 % diethylpyrocarbonate (DEPC) solution for 12 h at 37 °C before autoclaving or heating to 100 °C for 15 min to remove residual DEPC.
- Electrophoresis tanks should be cleaned with detergent solution (e.g., 0.5 % SDS), thoroughly rinsed with RNase-free water, rinsed with ethanol and finally allowed to dry.
- All buffers must be prepared with DEPC-treated RNase-free double-distilled water.
- Avoid handling bacterial cultures, cell cultures or other biological sources of RNases in the same lab where the RNA purification will be performed.
- Do not use equipment, glassware and plastic ware employed for other applications which might introduce RNase contaminations in the RNA isolation.


GenUP™ Total RNA Kit

SHORT PROTOCOL


STEPS

SCHEME


- Lyse the sample material.
- Centrifuge to pellet unlysed material.




- Transfer the lysate to the Mini-Filter (blue ring).
- Centrifuge, and discard the Mini-Filter with bound DNA.
- The filtrate contains the RNA, do not discard.




- Add ethanol to the filtrate.
- Transfer to a Mini-Filter (violet ring) to bind RNA and centrifuge.



- Add Buffer WASH A and centrifuge.
- Add Buffer WASH B and centrifuge.
- Centrifuge once more to remove residual ethanol.



- Add Water, RNase-free (for ELUTION), incubate and centrifuge.
- RNA in the Elution Tube is ready for use.



PROTOCOL FOR ISOLATION OF TOTAL RNA FROM TISSUE SAMPLES

PROCEDURE	NOTES
HOMOGENIZATION WITH A ROTOR-STATOR HOMOGENIZER - Transfer up to 20 mg fresh or frozen starting material to a suitable reaction vessel for the homogenizer. - Add 450 µl Buffer LYSIS LR and homogenize the sample. - Transfer the homogenized tissue sample to a 1.5 ml reaction tube.	Incomplete homogenization can reduce RNA yield.
HOMOGENIZATION WITH A MORTAR, PESTLE AND LIQUID NITROGEN - Transfer up to 20 mg fresh or frozen starting material to a mortar containing liquid nitrogen and grind to a fine powder. - Transfer the powder into a 1.5 ml reaction tube. Do not allow the sample to thaw. - Add 450 µl Buffer LYSIS LR and incubate the sample under continuous shaking until it is lysed completely (lysate becomes clear).	
Centrifuge at maximum speed for 1 min to pellet unlysed material.	
- Transfer the supernatant to a Mini-Filter (blue) placed in a Collection Tube. - Discard the reaction tube. - Centrifuge at $10,000 \times g$ (~12,000 rpm) for 2 min. - Discard the Mini-Filter with bound DNA and keep the filtrate.	- If the solution does not pass completely through the Mini-Filter, repeat the centrifugation or use a longer centrifugation time. - Do not discard the Collection Tube containing the RNA.
Add an equal volume of 70 % ethanol (400 µl) to the filtrate and mix by pipetting.	
- Transfer the sample to a Mini-Filter (violet) placed in a new Collection Tube. - Centrifuge at $10,000 \times g$ (~12,000 rpm) for 2 min. - Discard the Collection Tube with the filtrate.	If the solution does not pass completely through the Mini-Filter, repeat the centrifugation or use a longer centrifugation time.
- Place the Mini-Filter into a new Collection Tube. - Add 500 µl Buffer WASH A to the Mini-Filter. - Centrifuge at $10,000 \times g$ (12,000 rpm) for 1 min. - Discard the Collection Tube with the filtrate.	Before use, prepare Buffer WASH A as described above.

GenUP™ Total RNA Kit

- Place the Mini-Filter into a new Collection Tube.
 - Add 700 µl Buffer WASH B.
 - Centrifuge at 10,000 × g (~12,000 rpm) for 1 min.
 - Discard the Collection Tube with the filtrate.
-
- Place the Mini-Filter into a new Collection Tube.
 - Centrifuge at maximum speed for 2 min to remove residual ethanol.
 - Discard the Collection Tube.
-
- Place the Mini-Filter with bound RNA into an Elution Tube.
 - Add 30–80 µl Water, RNase-free (for ELUTION), to the center of the Mini-Filter RNA.
 - Incubate at room temperature for 1 min.
 - Centrifuge at 6,000 × g (8,000 rpm) for 1 min.
 - Discard the Mini-Filter.
-
- Purified RNA in the Elution Tube can be used immediately.
 - To improve yield, perform elution twice using ½ volume of Water, RNase-free (for ELUTION). Elute with at least 20 µl.
 - Store the RNA at 4 °C (short-term) or –80 °C (long-term).
-

PROTOCOL FOR ISOLATION OF TOTAL RNA FROM EUKARYOTIC CELLS

PROCEDURE	NOTES
• Transfer up to max 5×10^6 cells to an appropriate reaction tube and pellet by centrifugation. • Discard the supernatant. • Resuspend the cells in 400 μl Buffer LYSIS LR. • Incubate at room temperature for 2 min. • Resuspend by carefully pipetting up and down, and incubate at room temperature for an additional 3 min.	• Incomplete disruption can reduce RNA yield. • No cell clumps should be visible after lysis.
• Transfer the lysate to a Mini-Filter (blue) placed in a Collection Tube. • Discard the reaction tube. • Centrifuge at $10,000 \times g$ (~12,000 rpm) for 2 min. • Discard the Mini-Filter with bound DNA and keep the filtrate.	• If the solution does not pass completely through the Mini-Filter, repeat the centrifugation or use a longer centrifugation time. • Do not discard the Collection Tube containing the RNA.
• Add an equal volume of 70 % ethanol (400 μl) to the filtrate and mix by pipetting.	
• Transfer the sample to a Mini-Filter (violet) placed in a new Collection Tube. • Centrifuge at $10,000 \times g$ (~12,000 rpm) for 2 min. • Discard the Collection Tube with the filtrate.	• If the solution does not pass completely through the Mini-Filter, repeat the centrifugation or use a longer centrifugation time.
• Place the Mini-Filter with bound RNA into a new Collection Tube. • Add 500 μl Buffer WASH A to the Mini-Filter. • Centrifuge at $10,000 \times g$ (12,000 rpm) for 1 min. • Discard the Collection Tube with the filtrate.	• Before use, prepare Buffer WASH A as described above.
• Place the Mini-Filter into a new Collection Tube. • Add 700 μl Buffer WASH B. • Centrifuge at $10,000 \times g$ (~12,000 rpm) for 1 min. • Discard the Collection Tube with the filtrate.	• Before use, prepare Buffer WASH B as described above.
• Place the Mini-Filter into a new Collection Tube. • Centrifuge at maximum speed for 2 min to remove residual ethanol. • Discard the Collection Tube.	

GenUP™ Total RNA Kit

- Place the Mini-Filter with bound RNA into an Elution Tube.
- Add 30–80 µl Water, RNase-free (for ELUTION), to the center of the Mini-Filter RNA.
- Incubate at room temperature for 1 min.
- Centrifuge at 6,000 × g (8,000 rpm) for 1 min.
- Discard the Mini-Filter.
- To improve yield, perform elution twice using ½ volume of Water, RNase-free (for ELUTION). Elute with at least 20 µl.
- Purified RNA in the Elution Tube can be used immediately.
- Store the RNA at 4 °C (short-term) or –80 °C (long-term).

PROTOCOL FOR ISOLATION OF TOTAL RNA FROM BACTERIAL CELLS

PROCEDURE	NOTES
• Pellet up to 1×10⁹ bacterial cells by centrifugation at 5000 × g (6000 rpm) for 2–5 min. • Completely remove the supernatant completely, removing drops with a pipette if necessary. • Resuspend the cell pellet completely in 100 µl TE buffer by pipetting up and down. Avoid foaming. • Add 5–10 µl (Gram-positive) or 1–2 µl (Gram-negative) 50 mg/ml lysozyme to the cell suspension. • Pipette carefully up and down until the solution becomes clear.	• Incomplete homogenization can reduce RNA yield. • Before use, prepare the lysozyme and TE buffer as described above. • The optimal amount of lysozyme and incubation time varies depending on cell type.
• Add 450 µl Buffer LYSIS LR to the clarified sample. • Resuspend by carefully pipetting up and down. • Incubate at room temperature for an additional 3 min.	• Incomplete lysis can reduce RNA yield. After lysis, lysate should be clear or viscous, with no cell clumps.
• Transfer the lysate to a Mini-Filter (blue) placed in a Collection Tube. • Discard the reaction tube. • Centrifuge at 10,000 × g (~12,000 rpm) for 2 min. • Discard the Mini-Filter with bound DNA and keep the filtrate.	• If the solution does not pass completely through the Mini-Filter, repeat the centrifugation or use a longer centrifugation time. • Do not discard the Collection Tube containing the RNA.
• Add an equal volume of 70 % ethanol (400 µl) to the filtrate and mix by pipetting.	
• Transfer the sample to a Mini-Filter (violet) placed in a new Collection Tube. • Centrifuge at 10,000 × g (~12,000 rpm) for 2 min. • Discard the Collection Tube with the filtrate.	• If the solution does not pass completely through the Mini-Filter, repeat the centrifugation or use a longer centrifugation time.

-
- | | |
|---|---|
| <ul style="list-style-type: none"> • Place the Mini-Filter RNA into a new Collection Tube. • Add 500 µl Buffer WASH A to the Mini-Filter. • Centrifuge at 10,000 × g (12,000 rpm) for 1 min. • Discard the Collection Tube with the filtrate. | <ul style="list-style-type: none"> • Before use, prepare Buffer WASH A as described above. |
|---|---|
-
- | | |
|---|---|
| <ul style="list-style-type: none"> • Place the Mini-Filter into a new Collection Tube. • Add 700 µl Buffer WASH B. • Centrifuge at 10,000 × g (~12,000 rpm) for 1 min. • Discard the Collection Tube with the filtrate. | <ul style="list-style-type: none"> • Before use, prepare Buffer WASH B as described above. |
|---|---|
-
- | | |
|--|--|
| <ul style="list-style-type: none"> • Place the Mini-Filter RNA into a new Collection Tube. • Centrifuge at maximum speed for 2 min to remove residual ethanol. • Discard the Collection Tube. | |
|--|--|
-
- | | |
|--|--|
| <ul style="list-style-type: none"> • Place the Mini-Filter with bound RNA into an Elution Tube. • Add 30–80 µl Water, RNase-free (for ELUTION), to the center of the Mini-Filter RNA. • Incubate at room temperature for 1 min. • Centrifuge at 6,000 × g (8,000 rpm) for 1 min. • Discard the Mini-Filter. | <ul style="list-style-type: none"> • To improve yield, perform elution twice using ½ volume of Water, RNase-free (for ELUTION). The minimum elution volume should exceed 20 µl. |
|--|--|
-
- | | |
|---|---|
| <ul style="list-style-type: none"> • Purified RNA in the Elution Tube can be used immediately. | <ul style="list-style-type: none"> • Store the RNA at 4 °C (short-term) or –80 °C (long-term). |
|---|---|
-

GenUP™ Total RNA Kit

TROUBLESHOOTING

PROBLEM	SOLUTION
CLOGGED MINI-FILTER	
Insufficient disruption or homogenization	Reduce amount of starting material. After lysis, centrifuge the lysate to pellet debris and continue with the protocol using the supernatant.
LOW YIELD	
Insufficient disruption or homogenization	Reduce amount of starting material. Avoid overloading the Mini-Filter, as overloading reduces yield.
Incomplete elution	To improve elution, prolong the incubation time to 5 min or repeat elution.
DNA CONTAMINATION	
Too much starting material	Reduce amount of starting material.
Incorrect lysis of starting material	Use the recommended techniques for lysis. Perform an on-column DNase digestion step after binding the RNA to the Mini-Filter RNA (violet). Alternatively, perform DNase digestion of the eluate. Ensure that the DNase I is RNase-free.
DEGRADED RNA	
RNA source inappropriately handled or stored	Ensure that the starting material is fresh. Ensure that the protocol – especially the first steps – has been performed quickly.
RNase contamination of solutions, tubes, etc.	Use sterile, RNase-free filter tips. Before every RNA preparation, clean the pipette, devices and working place. Always wear gloves.
RNA DOES NOT PERFORM WELL IN DOWNSTREAM APPLICATIONS (E.G., RT-PCR)	
Ethanol carryover	Increase centrifugation time for removing ethanol.
Salt carryover	Ensure that Buffer WASH A and Buffer WASH B are at room temperature before use. If a buffer contains precipitate, dissolve the precipitate by warming carefully.

SAFETY PRECAUTIONS

- This kit is made for single use only!
- Don't eat or drink components of the kit!
- The kit shall only be handled by educated personnel in a laboratory environment!
- Wear gloves while handling these reagents and avoid skin contact! In case of contact, flush with water immediately!
- Handle and discard waste according to local safety regulations!
- Do not add bleach or acidic components to the waste after sample preparation!
- Buffers LYSIS LR and WASH A contain guanidine isothiocyanate, which is harmful to health. Contact with acids liberates very toxic gas!

CERTIFICATE OF ANALYSIS

The kit was tested by isolation of RNA from tissue samples and subsequent analysis of RNA.

Quality confirmed by: Head of Quality Control

SAFETY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS):

Sicherheitshinweise finden Sie in den Sicherheitsdatenblättern (SDB) unter:

<http://www.biotechrabbit.com/support/documentation.html>

USEFUL HINTS

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

CONTACT BIOTECHRABBIT

biotechrabbit GmbH
Volmerstr. 9
12489 Berlin, Germany

info@biotechrabbit.com
support@biotechrabbit.com
www.biotechrabbit.com

Phone: +49 30 555 7821 10
Fax: +49 30 555 7821 99



LEGAL DISCLAIMER

Purchase of product does not include a license to perform any patented or otherwise protected applications; therefore, it is the sole responsibility of customer to determine whether any licenses for the particular application in which the product is used are required and to engage in the relevant license agreement(s).

LIMITATION OF USE

For Research Use Only. biotechrabbit does not assume any warranty regarding the design, merchantability or fitness for a particular purpose of the product, especially, but not limited to, for any *in vitro* diagnostic or other usages of the products. It is the sole responsibility of the customer to obtain the required approvals from the competent authorities. The product is not suitable for administration to humans or animals.

Trademarks: *biotechrabbit™*, *GenUP™* (biotechrabbit GmbH)

valid from 31.10.2025