

# GenUP™ Plant RNA Kit

LOT: See product label

**EXPIRY DATE:** See product label

# **ORDERING INFORMATION**

PRODUCT	GenUP™ Plant RNA Kit		
CAT.NO.	BR0701501	BR0701502	BR0701503
SIZE	10 preps	50 preps	250 preps
COMPONENTS			
Buffer LYSIS LR	15 ml	30 ml	160 ml
Buffer LYSIS LT	8 ml	30 ml	130 ml
Buffer WASH A (concentrate)	5 ml (add 5 ml ethanol)	15 ml (add 15 ml ethanol)	70 ml (add 70 ml ethanol)
Buffer WASH B (concentrate)	6 ml (add 24 ml ethanol)	16 ml (add 64 ml ethanol)	2 × 36 ml (add 144 ml ethanol)
Water,RNase-free (for ELUTION)	2ml	3×2 ml	25 ml
Mini Filters DNA (blue)	10	50	5×50
Mini Filters RNA (violet)	10	50	5×50
Collection Tubes (2.0 ml)	60	6×50	30×50
Elution Tubes (1.5 ml)	10	50	5×50

**STORAGE** 

Room temperature (until expiry date - see product label).

If precipitation appears, gently warm the solution to dissolve the precipitate.

# **FEATURES**

- Fast and simple procedure
- High-quality RNA isolated from a wide variety of plant samples
- Physical removal of DNA, no DNase treatment, no toxic β-mercaptoethanol

# **APPLICATIONS**

• Isolation of total RNA from up to 100 mg plant material

### DESCRIPTION

biotechrabbit™ GenUP Plant RNA Kit has been developed for quick and easy purification of total RNA from plant materials. After initial homogenization and lysis, genomic DNA is bound to a Mini Filter DNA, which can be discarded. RNA is selectively bound to a Mini Filter RNA, washed with two different buffers and eluted. The purified RNA is ready for use in any demanding molecular biology application, including RT-PCR. Two lysis buffers, Buffer LYSIS LR and Buffer LYSIS LT, are provided to maximize yield. Most plant material can be processed with Buffer LYSIS LR. In the cases that yield using Buffer LYSIS LR is low, use Buffer LYSIS LT.

# **SPECIFICATIONS**

STARTING MATERIAL	Plant material (up to 100 mg)
EXTRACTION TIME	30 min after homogenization
BINDING CAPACITY	Approximately 100 µg RNA
TYPICAL YIELD	Variable depending on the starting material; approximately 70 $\mu gRNA$
AVERAGE PURITY	A260/A280 1.7-2.0

## MATERIALS SUPPLIED BY THE USER

- 96–99.8% ethanol
- 70% ethanol
- Reaction tubes
- Pipet tips
- Optional: DNase I

## STEPS BEFORE STARTING

Add the following volume of 96–99.8% ethanol to the indicated buffer bottle, close firmly, mix thoroughly
and store at room temperature.

CAT.NO.		CONCENTRATE	ETHANOL	FINAL VOLUME
Buffer WASH A	BR0701501	5 ml	5ml	10 ml
	BR0701502	15 ml	15 ml	30 ml
	BR0701503	70 ml	70 ml	140 ml
Buffer WASH B	BR0701501	6 ml	24 ml	30 ml
	BR0701502	16 ml	64 ml	80 ml
	BR0701503	36 ml	144 ml	180 ml

- Avoid repeated freezing and thawing of starting materials.
- Mark all vials and filters to avoid confusion when purifying multiple preps.
- Centrifugation steps should be carried out at room temperature.

### GUIDELINES FOR PREVENTION OF RNA DEGRADATION


Special care should be taken to minimize contamination with RNases, as RNA is extremely sensitive to degradation.

- · Always wear gloves, and change them frequently.
- Keep all tubes closed when possible.
- Keep samples and isolated RNA on ice.
- Reduce preparation time as much as possible.
- Use only sterile, disposable polypropylene tubes throughout the procedure (these tubes are generally RNase-free).
- Non-disposable plastic ware should be treated before use to ensure that it is RNase-free. Plastic ware should be thoroughly rinsed with 0.1 M NaOH, 1 mM EDTA followed by RNase-free water. You can also take chloroform-resistant plastic ware rinsed with chloroform to inactivate RNases.
- All glassware should be treated before use to ensure that it is RNase-free.
  - Glassware should be cleaned with detergent, thoroughly rinsed and oven baked at 240 °C for four or more hours before use. Oven baking inactivates RNases and ensures that no other nucleic acids (such as plasmid DNA) are present on the surface of the glassware.
  - Autoclaving alone will not inactivate many RNases completely. The glassware should be immersed in 0.1% diethylpyrocarbonate (DEPC) solution for 12 h at 37°C before autoclaving or heating to 100°C for 15 min to remove residual DEPC.
- Electrophoresis tanks should be cleaned with detergent solution (e.g., 0.5% SDS), thoroughly rinsed with RNase-free water, rinsed with ethanol and finally allowed to dry.
- All buffers must be prepared with DEPC-treated RNase-free double-distilled water.
- Avoid handling bacterial cultures, cell cultures or other biological sources of RNases in the same lab where
  the RNA purification will be performed.
- Do not use equipment, glassware and plastic ware employed for other applications which might introduce RNase contaminations in the RNA isolation.


# SHORT PROTOCOL

STEPS SCHEME


- Homogenize plant material and lyse with Buffer LYSIS LR or Buffer LYSIS LT
- Centrifuge to remove unlysed material.




- Transfer the lysate to the Mini Filter DNA (blue).
- · Centrifuge, and discard the Mini Filter DNA.
- The filtrate contains the RNA, do not discard.




- Add ethanol to the filtrate.
- Transfer to a Mini Filter RNA (violet) to bind RNA and centrifuge.



- Add Buffer WASH A and centrifuge.
- Add Buffer WASH B and centrifuge.
- Centrifuge once more to remove residual ethanol.



- Elute RNA with Water, RNase-free (for ELUTION), and centrifuge.
- Purified RNA in the Elution Tube is ready for use.



# PROTOCOL FOR ISOLATION OF TOTAL PLANT RNA

#### **PROCEDURE NOTES** • Homogenize up to 100 mg plant material completely. • For homogenization, use liquid nitrogen with a pestle and mortar or a rotor-stator • Transfer 450 ul Buffer LYSIS LR or Buffer LYSIS LT to a homogenizer. When using liquid nitrogen, do 1.5 ml reaction tube. not let the sample thaw. When using a rotor-• Transfer the homogenized material to the tube containing stator homogenizer, the plant material can be the lysis buffer. homogenized in the lysis buffer. The lysate can be stored at -20°C. Centrifuge at maximum speed for 1 min to pellet debris. • Transfer the supernatant to a Mini Filter DNA (blue) • If the solution does not pass completely placed in a Collection Tube. through the Mini Filter, repeat the Discard the reaction tube. centrifugation or use a longer centrifugation Centrifuge at 10.000 × g (12.000 rpm) for 2 min. Do not discard the Collection Tube containing Discard the Mini Filter DNA. the RNA Add an equal volume of 70% ethanol (400 µl) to the filtrate and mix by pipetting. • Transfer sample to the Mini Filter RNA (violet) placed in a If the solution does not pass completely new Collection Tube. through the Mini Filter, repeat the centrifugation or use a longer centrifugation Centrifuge at 10,000 × g (12,000 rpm) for 2 min. time. Discard the Collection Tube with the filtrate. Place the Mini Filter RNA into a new Collection Tube. • Before use, prepare Buffer WASH A as Add 500 ul Buffer WASH A to the Mini Filter. described above. Centrifuge at 10.000 × q (12.000 rpm) for 1 min. Discard the Collection Tube with the filtrate. Place the Mini Filter BNA into a new Collection Tube. Before use, prepare Buffer WASH Bas described above. Add 650 µl Buffer WASH B to the Mini Filter. Centrifuge at 10,000 × g (12,000 rpm) for 1 min. Discard the Collection Tube with the filtrate. Place the Mini Filter RNA into a new Collection Tube. Add 650 ul Buffer WASH B to the Mini Filter. • Centrifuge at 10,000 × g (12,000 rpm) for 1 min. Discard the Collection Tube with the filtrate. Place the Mini Filter RNA into a new Collection Tube. • Centrifuge at maximum speed for 2 min to remove residual ethanol. · Discard the Collection Tube. Place the Mini Filter BNA into a new Flution Tube. • To improve yield, perform elution twice using • Add 30-80 µl Water, RNase-free (for ELUTION). 1/2 volume of Water, RNase-free (for ELUTION). The minimum elution volume Incubate at room temperature for 1 min. should exceed 20 µl. • Centrifuge at 6000 × g (8000 rpm) for 1 min. Discard the Mini Filter RNA.

Purified RNA in the Elution Tube can be used immediately.

Store the RNA at 4°C (short-term) or -20°C

(long-term).

# GenUP™ Plant RNA Kit

TROUBLESHOOTING			
PROBLEM	SOLUTION		
CLOGGED MINI FILTER			
Insufficient disruption or homogenization	Reduce the amount of starting material.  After lysis, centrifuge lysate to pellet debris and continue the protocol using the supernatant.		
LOWYIELD			
Insufficient disruption or homogenization	Reduce the amount of starting material. Do not overload the Mini Filter DNA.		
Incomplete elution	Increase the elution time up to 5 min, or repeat the elution. Use a higher elution volume or elute in two steps.		
DNA CONTAMINATION			
Too much starting material	Reduce the amount of starting material.		
Incorrectlysis of starting material	Use the recommended technique to lysis the starting material.  Perform an on-filter DNase digestion after binding RNA to the Mini Filter RNA.  Alternatively, perform DNase digestion on the eluate.		
TOTAL RNA DEGRADED			
Starting material handled or stored in appropriately	Ensure the starting material is fresh. Perform the homogenization and lysis steps quickly.		
RNase contamination	Use sterile, RNase-free filter tips. Ensure pipets, devices and workplace are clean. Wear gloves.		
PURIFIED RNA DOES NOT PERFORM WELL IN DO WNSTREAM APPLICATIONS (E.G., RT-PCR)			
Ethanol carried over to elution	Increase the centrifugation time when removing ethanol.		
Salt carried over to elution	Ensure that Buffer WASH A and Buffer WASH B are at room temperature. Dissolve precipitates in the washing solutions by warming carefully.		

### SAFETY PRECAUTIONS

- This kit is made for single use only!
- Don't eat or drink components of the kit!
- The kit shall only be handled by educated personal in a laboratory environment!
- Wear gloves while handling these reagents, and avoid skin contact! In case of contact, flush with water immediately!
- Handle and discard waste according to local safety regulations!
- Do not add bleach or acidic components to the waste after sample preparation!
- Buffers LYSIS LR and WASH A contain guanidine isothiocyanate, which is harmful to health. Contact with acids liberates very toxic gas!

### **CERTIFICATE OF ANALYSIS**

The components of the kit were tested for isolation of total RNA from plant material and subsequent analysis.

Quality confirmed by: Head of Quality Control

# **SAFETY INSTRUCTIONS**

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter: http://www.biotechrabbit.com/support/documentation.html.

# **USEFUL HINTS**

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

## CONTACT BIOTECHRABBIT

biotechrabbit GmbH info@biotechrabbit.com

Volmerstr. 9a support@biotechrabbit.com Phone: +49 30 555 7821-10

12489 Berlin. Germany www.biotechrabbit.com Fax: +49 30 555 7821-99



### Legal Disclaimer and Product Use Limitation

Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.

Trademarks: biotechrabbit™, GenUP™ (biotechrabbit GmbH).

valid from 24.08.2016