

GenUP™ FFPE Paraffin Removal Solution

LOT: See product label EXPIRY DATE: See product label

ORDERING INFORMATION

CAT. NO.	SIZE	PACKAGE CONTENT
BR0701601	10 ml	10 ml GenUP FFPE Paraffin Removal Solution
BR0701602	50 ml	5 × 10 ml GenUP FFPE Paraffin Removal Solution

COMPONENT	COMPOSITION
GenUP FFPE Paraffin Removal Solution	Single solution formulation in dispensing-cap bottle for fast and easy DNA or RNA extraction.
STORAGE	Room temperature (until expiry date – see product label)

FEATURES

- One-step deparaffinization of FFPE slices
- No toxic or flammable chemicals

APPLICATIONS

- Removal of paraffin from formalin-fixed, paraffin-embedded (FFPE) sample slices
- Deparaffinization for downstream DNA and RNA purification

SPECIFICATIONS

STARTING MATERIAL	2–4 fresh slices from a FFPE sample
EXTRACTION TIME	Approximately 45–60 min
YIELD	Depends on the type, quality and amount of starting material Typically, 50-200 reactions per 10 ml FFPE
DNA AND RNA ISOLATION	Compatible with all common silica-based DNA and RNA purification kits

GenUP™ FFPE Paraffin Removal Solution

DESCRIPTION

biotechrabbit™ GenUP FFPE Paraffin Removal Solution is an odor-free, biofriendly and innovative solution to deparaffinize FFPE slices within minutes in a single step. This product replaces conventional procedures that involve toxic and/or flammable solvents. After removal of paraffin from FFPE slices with GenUP FFPE Paraffin Removal Solution, continue with the genomic DNA or RNA purification kit of your choice.

MATERIALS SUPPLIED BY THE USER

- · Reaction tubes
- Pipet tips
- GenUP gDNA Kit (BR0700601) or other nucleic acid purification kit
- GenUP™ Total RNA Kit (BR0700901) or other nucleic acid purification kit

STEPS BEFORE STARTING

- The GenUP FFPE Paraffin Removal Solution has a freezing point of 18°C. Buffer in a solid aggregate state can be heated to 25–35°C in a water bath to melt.
- Heat a thermomixer to 60–65°C.
- Place the buffers, tubes and columns from the nucleic acid purification kit on the bench so they will be available for immediate use after the paraffin removal step.


SHORT PROTOCOL

STEPS SCHEME

 Cover FFPE sample slices with GenUP FFPE Paraffin Removal Solution, add the lysis buffer from your genomic DNA purification or total RNA extraction kit and vortex.



 Incubate first at 60–65°C to degrade proteins. Then incubate for an additional 10 min at 90°C.



 Transfer the lower phase to a new tube and continue with the GenUP gDNA, GenUP Total RNA or other nucleic acid purification kit.



PROTOCOL

PROCEDURE	NOTES
Cut FFPE slices from your FFPE sample with a microtome.	
Place 2–4 FFPE slices in a 1.5 ml reaction tube.	Depending on the tissue and the age of the FFPE block, more slices are needed in the case of low recovery.
 Add 3–5 drops of GenUP FFPE Paraffin Removal Solution to the FFPE slices. 	• The FFPE slices dissolve immediately. Add enough solution to submerge the FFPE slices (approximately 300–500 μl).
Add buffers for the first step of GenUP DNA or RNA purification kit.	
Vortex the reaction mixture for 10 s.	
 Place the reaction mixture in the thermomixer and incubate at 60–65°C for the time indicated in the purification kit protocol. 	 Incubate at least for 15–30 min until FFPE slices are dissolved.
 Incubate for an additional 10 min at 90°C. Remove the reaction mixture from the thermomixer. 	If the two phases have not separated, centrifuge the reaction mixture briefly to collect the aqueous phase at the bottom of the reaction tube.
Transfer the lower colorless aqueous phase containing the DNA or RNA to a new 1.5 ml reaction tube.	The upper, orange phase contains the paraffin.

 Continue with the protocol for the GenUP gDNA or Total RNA Kit. For RNA purification, start with adding ethanol.

SAFETY PRECAUTIONS

- This kit is made for single use only!
- Don't eat or drink components of the kit!
- The kit shall only be handled by educated personal in a laboratory environment!
- Wear gloves while handling these reagents, and avoid skin contact! In case of contact, flush with water immediately.
- · Handle and discard waste according to local safety regulations.
- Do not add bleach or acidic components to the waste after sample preparation!

GenUP™ FFPE Paraffin Removal Solution

CERTIFICATE OF ANALYSIS

A distinct band of 750 bp was amplified using specific primers after DNA extraction from a paraffin block with embedded human tissue by PCR.

Quality confirmed by: Head of Quality Control

SAFETY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter: http://www.biotechrabbit.com/support/documentation.html.

USEFUL HINTS

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

CONTACT BIOTECHRABBIT

biotechrabbit GmbH

Volmerstr. 9a info@biotechrabbit.com
12489 Berlin, support@biotechrabbit.com Phone: +49 30 555 7821-10
Germany www.biotechrabbit.com Fax: +49 30 555 7821-99



Legal Disclaimer and Product Use Limitation

Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.

Trademarks: biotechrabbit™, GenUP™ (biotechrabbit GmbH). valid from 30.08.2019