

Long Range DNA Polymerase, 2.5 U/ μ l

LOT: See product label

EXPIRY DATE: See product label

ORDERING INFORMATION

CAT. NO.	SIZE	PACKAGE CONTENT
BR0300301	100 U (80 rxn)	40 μ l Long Range PCR Enzyme Mix 400 μ l 10 \times Long Range Reaction Buffer 1.5 ml 5 \times PCR Enhancer
BR0300302	500 U (400 rxn)	5 \times 40 μ l Long Range PCR Enzyme Mix 5 \times 400 μ l 10 \times Long Range Reaction Buffer 1.5 ml 5 \times PCR Enhancer*

COMPONENT

COMPOSITION

Long Range PCR Enzyme Mix	Long Range PCR Enzyme Mix, 2.5 U/ μ l, in storage buffer containing 50% (v/v) glycerol
10 \times Long Range Reaction Buffer	Optimized PCR reaction buffer including magnesium ions
5 \times PCR Enhancer*	Proprietary PCR enhancer mix

* For reaction optimization. Sufficient PCR Enhancer for all reactions to be ordered separately (BR1900201).

STORAGE

-20°C (until expiry date – see product label)

FEATURES

- High-productivity, long-range PCR
- Increased fidelity for accurate amplification of GC-rich templates
- Polymerase mix for high yield and short cycle times

APPLICATIONS

- Long-range PCR up to 40 kb
- Amplification of GC-rich templates

Long Range DNA Polymerase, 2.5 U/ μ l

DESCRIPTION

biotechrabbit™ Long Range DNA Polymerase is a first-choice for amplification of targets up to 40 kb in size with higher accuracy than *Taq* DNA polymerase.

This specially designed blend of thermophilic polymerases is well suited for amplification of targets that are GC-rich and have complex structures.

For the most demanding applications, the supplied 5 \times PCR Enhancer can be optionally used to improve results when using templates with GC-rich sequences and complex structures.

Long Range DNA Polymerase produces a mixture of A-tailed and blunt-end PCR products. It is advisable to blunt products before cloning into blunt-end vectors.

PROTOCOL

Prevention of PCR contamination

When assembling the amplification reactions, care should be taken to eliminate the possibility of contamination with undesired DNA.

- Use separate clean areas for preparation of samples and reaction mixtures and for cycling.
- Wear fresh gloves. Use sterile tubes and pipette tips with aerosol filters for PCR setup.
- Use only water and reagents that are free of DNA and nucleases.
- With every PCR setup, perform a contamination control reaction that does not include template DNA.

Standard PCR setup

The standard PCR protocol using biotechrabbit reaction buffer provides excellent results for most applications. Optimization might be necessary for certain conditions, such as the amplification of long targets, high GC or AT content, strong template secondary structures or insufficient template purity. In such cases, optimization of template purification (see biotechrabbit nucleic acid purification kits), primer design and annealing temperature is recommended.

The best conditions for each primer-template can be optimized with the following:

- Choosing the optimal quantities of template and primers
- Determining optimal concentrations of the enzyme and magnesium ions
- Optimizing cycling conditions

If unspecific amplification occurs, the amount of DNA Polymerase and the primer concentration can be reduced. Correspondingly, these can be increased when yield is low.

BASIC PROTOCOL

- Optionally, use the supplied 5 \times PCR Enhancer to increase the yield and to lower the background in more complicated PCR reactions (low amounts of template, impure or GC-rich template).
- Thaw on ice and mix all reagents well, especially the $MgCl_2$ solution and dNTPs.
- Keep all reagents and reactions on ice.
- When setting up multiple reactions, prepare a master mix of water, buffer, dNTPs and polymerase. Prepare enough master mix for one more than the actual number reactions. Alternatively, use biotechrabbit Long Range PCR Master Mix, 2 \times (cat. no. BR0300401)
- Pipet the master mix into thin-walled 0.2 ml PCR tubes.
- Add template and primers separately if they are not used in all reactions.

Long Range DNA Polymerase, 2.5 U/ μ l

COMPONENT	VOLUME	FINAL CONCENTRATION
10 \times Long Range Reaction Buffer	5 μ l	1 \times
5 \times PCR Enhancer (optional)*	10 μ l	1 \times
* For reaction optimization. Sufficient PCR Enhancer for all reactions to be ordered separately (BR1900201).		
10 mM dNTP Mix	1 μ l	200 μ M
Forward primer	Variable	0.2–1 μ M
Reverse primer	Variable	0.2– μ M
Template DNA	Variable	10 pg–1 μ g
<i>Use 0.01–1 ng for plasmid or phage DNA and 0.1–1 μg for genomic DNA</i>		
Long Range PCR Enzyme Mix (2.5 U/ μ l)	0.5–1 μ l	1.25–2.5 U
Nuclease free water	Variable	
Total volume	50 μ l	

- Mix and centrifuge briefly to collect the liquid in the bottom of the tube.
- Place in the PCR cyclor.

CYCLING PROGRAM

STEP	TEMPERATURE	TIME	CYCLES
Initial activation	95°C	2 min	1
Denaturation	95°C	30 s	25–35
Annealing	55°C	30–45 s	25–35
<i>Approximately 5°C below T_m of primers</i>			
Extension	72°C	30 s/kb	25–35
<i>For fragments longer than 5 kb, use 68°C extension temperature and 1 min/kb timing.</i>			
Final extension	72°C	5 min	1
<i>To extend all incomplete PCR products</i>			
Storage in the cyclor	4°C	Indefinitely	1

- Add loading dye solution (see DNA Loading Dye, 6 \times , cat. no. BR0800301) to the reactions to analyze PCR products on a gel or store them at –20°C.
- For cloning, always purify the PCR product from a gel (see BR0700401 GenUP™ Gel Extraction Kit).

Long Range DNA Polymerase, 2.5 U/ μ l

CERTIFICATE OF ANALYSIS

Unit Definition

One unit catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 minutes at 72°C.

Quality Control

Functional assay

Human genomic DNA was amplified using the DNA Polymerase and specific primers to produce a distinct band of 750 bp.

Self-priming activity

Standard PCR is carried out without primers, using the DNA Polymerase and human genomic DNA. No products were amplified.

Exonuclease assay

Linearized lambda/HindII fragments are incubated with the DNA Polymerase in a 50 μ l reaction mixture for 4 h at 37°C. No degradation of DNA was observed.

Endonuclease assay

lambda DNA is incubated with the DNA Polymerase in a 50 μ l reaction mixture for 4 h at 37°C. No degradation of DNA was observed.

Nick Activity

Supercoiled plasmid DNA is incubated with the DNA Polymerase in a 50 μ l reaction mixture for 4 h at 37°C. No conversion of covalently closed circular DNA to nicked DNA was detected.

E. coli DNA contamination assay

A sample of the denatured DNA Polymerase is analyzed with specific primers targeting the 16S rRNA gene in qPCR for the presence of contaminating *E. coli* DNA. No *E. coli* DNA was detectable.

Quality confirmed by: Head of Quality Control

SAFETY INSTRUCTIONS

For safety instructions please see Safety Data Sheets (SDS)/Sicherheitshinweise finden Sie in den SDS unter: <http://www.biotechrabbit.com/support/documentation.html>.

USEFUL HINTS

- Visit Applications at www.biotechrabbit.com for more products and product selection guides.
- Most biotechrabbit products are available in custom formulations and bulk amounts.

CONTACT BIOTECHRABBIT

biotechrabbit GmbH
Volmerstr. 9a
12489 Berlin, Germany

info@biotechrabbit.com
support@biotechrabbit.com
www.biotechrabbit.com

Phone: +49 30 555 7821-10
Fax: +49 30 555 7821-99



Legal Disclaimer and Product Use Limitation

Purchase of product does not include a license to perform any patented applications; therefore it is the sole responsibility of users to determine whether they may be required to engage a license agreement depending upon the particular application in which the product is used. This product was developed, manufactured, and sold for in vitro use only. It is not suitable for administration to humans or animals.

Trademarks: biotechrabbit™ (biotechrabbit GmbH).

valid from 24.01.2017